NOVA Cer[®] LDS

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*Instruction of Use
Lithium Disilicate layering ceramic*

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Safety information: *Please read these instructions thoroughly before using the product. Improper use of the product and failure to observe the information provided may affect quality and reduce the lifetime of the dental prosthesis. Only trained staff (dentists or dental technicians) should be allowed to handle this product.*

– WARNING –

Applies to porcelains in powder form: Please open the packaging and use the product only with protective equipment. Wear safety goggles and a dust mask. Dust particles irritate the eyes, mucous membranes and skin, and injure the lungs.

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Intended use: NOVACer® LDS is a veneering ceramic for the veneering of Lithium Disilicate frameworks. When veneering Lithium Disilicate frameworks, the guidelines of the respective manufacturer must be observed.

The framework should have at least 50% of the total thickness.

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BUILD UP PROCEDURE:

Frame:

The lithium disilicate copings were pressed or produced in the CAD/CAM procedure.

Note: *For lithium disilicate framework dimension manufacturer's instruction must be respected!*

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The labial and occlusal shape is build-up with NOVACer® LDS dentine in the required shade.

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*The dentin will be reduced and chromatic mamelon
material build-up.*

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*Control of the correct volume of the
dentin structure after firing.*

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*The missing areas will be build-up with
NOVACer[®] LDS incisal and tranparent porcelains.*

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Glaze firing after finishing.

The degree of glaze can be adjusted by polishing.

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* The firing parameters given above are guidelines, which must always be adjusted to suit the furnace used and the situation of the furnace. What is essential is getting the right firing result.

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*All tested materials conform to EN ISO 9693-2. The technical and
physical values quoted relate to samples produced in-house at our
own measuring instruments.*

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